

Report No.: S-2016-03411 AM
Version: English
Page: 1 of 39
Print date: September 24th, 2016

Final Report S-2016-03411 AM

VALIDATION OF A GC/MS METHOD FOR THE IDENTIFICATION AND QUANTIFICATION OF NICOTINE RESIDUE IN THE TEST ITEM "VX1" AND ANALYSIS ON FIVE PRODUCTION BATCHES

Study program: S-2016-03411 AM

Contract n: QCH1071

Sponsor: [REDACTED]

Test facility: [REDACTED]

Test item: "VX1"

Study Director: Released on: Sep 24th, 2016

INDEX

INDEX	2
COMPLIANCE WITH GOOD LABORATORY PRACTICE	4
QUALITY ASSURANCE STATEMENT	5
SUMMARY	6
INTRODUCTION	6
BIBLIOGRAPHY	6
FILING	7
PROCEDURES	7
TEST ITEM IDENTIFICATION	7
REFERENCE STANDARD	9
EXPERIMENTATION S-2016-03411 AM: GC/MS METHOD FOR THE QUANTIFICATION OF NICOTINE AND FIVE BATCHES ANALYSIS	10
INFORMATION	10
EXPERIMENTAL PROCEDURE - VALIDATION	10
TEST METHOD	10
<i>Parameters under investigation</i>	10
<i>Acceptability criteria</i>	11
<i>Analytical sequence</i>	12
<i>Analytical acceptability criteria</i>	12
<i>Equipment</i>	12
<i>Reagent</i>	13
<i>Working solutions preparation</i>	13
GC-MS METHOD – SIM MODE	15
CALCULATIONS	16
EXPERIMENTAL PROCEDURE - GC/MS CONFIRMATORY TECHNIQUE	17
TEST METHOD – GC/MS	17
<i>Equipment</i>	17
<i>Analysis</i>	17
<i>Working solutions preparation</i>	17
RESULTS	19
VALIDATION METHOD PARAMETERS	19
<i>Selectivity</i>	19
<i>Specificity</i>	21
<i>Linearity</i>	23
<i>Accuracy</i>	24
<i>Precision</i>	24
<i>Repeatability</i>	25
<i>LOQ</i>	25
<i>ANOMALY TEST (DIXON, GRUBBS)</i> :	25
<i>Statistical consideration</i>	25
DETERMINATION OF (-)-NICOTINE CONTENT IN FIVE BATCHES OF THE TEST ITEM "VX1"	26
DEVIATION	26
CONCLUSIONS	26
ANNEXES	26
ANNEX#1: NICOTINE - REFERENCE STANDARD CoA	27
ANNEX#2 : LINEARITY - EXCEL SHEET	30
ANNEX#3 : PRECISION-REPEATABILITY - EXCEL SHEET	31
ANNEX#4 : ACCURACY AND 5 BATCHES ANALYSIS - EXCEL SHEET	33
ANNEX#5 : STATISTIC - EXCEL SHEET	35

		Report No.: S-2016-03411 AM
		Version: English
		Page: 3 of 39
		Print date: Sept 24 th , 2016

ANNEX#6 : LOQ - EXCEL SHEET	36
ANNEX#7 : ANOMALY TEST - EXCEL SHEET	37
1) <i>Shapiro-Wilk Test</i>	37
2) <i>Grubbs Test</i>	38
3) <i>Dixon Test</i>	39

COMPLIANCE WITH GOOD LABORATORY PRACTICE

I the undersigned declare that the studies described in this report have been conducted under my supervision and in compliance with the following standards of Good Laboratory Practice:

- OECD Series on Principles of Good Laboratory Practice and Compliance Monitoring - OECD principles of Good Laboratory Practice (as revised in 1997) – Environment Directorate – Organisation for Economic Co- Operation and Development, Paris 1998.
- Legislative decree n. 50 of March the 2nd, 2007. Enforcement of Community Directives 2004/9/CE e 2004/10/CE, concerning the inspection and verification of Good Laboratory Practice and the drawing of the legislative, regulatory and administrative dispositions relative to the application of Good Laboratory Practice rules, to the control of their application on the assays performed on the chemical substances (GU n.86 of April the 13th, 2007).
- United States Food and Drug Administration, Title 21 Code of Federal Regulations Part 58, Federal Register 22 December 1978, and subsequent amendments.
- Certification N. 038/2013 released by the Italian Ministry of Health on November 19th 2013 and Provisional Certificate released on November 20th 2015 authorizing [REDACTED] to perform analyses in compliance with the principles of good laboratory practices

There were no circumstances that may affected the quality or integrity of the study.

[REDACTED]
Study Director

Sept 24th 2016
Date

QUALITY ASSURANCE STATEMENT

The study was assessed for compliance with the approved study program and the Standard Operating Procedures of [REDACTED]

The study and/or the test facility were periodically inspected by the Quality Assurance unit according to the corresponding SOPs. These inspections and audit were carried out by the Quality Assurance unit, personnel independent of staff involved in the study.

The undersigned hereby certifies the dates on which the inspections have been carried out and reported to the Director of the Study and to [REDACTED] Management:

QAU INSPECTIONS	
PHASE	DATE
Experimentation: -Audit process-based <i>Validation of analytical methods (GC)</i>	October 20 th - 23 th 2015
<i>Assay determination according to validated method (GC)</i>	June, 07 th 2016
-Audit study-based	//
Documentation: - Study program - Raw data - Final report	September, 24 th 2016 September, 24 th 2016 September, 24 th 2016

This report accurately reflects the raw data.

[REDACTED]
QA GLP
[REDACTED]

Sept 24th 2016
DATE

SUMMARY

The aim of the study was to set up and validate a GC/MS method for the identification and the quantification of the impurity (-)-Nicotine in the test item "VX1 validation" and then test this analyte in 5 different production batches of the test item "VX1".

The GC/MS-SIM Mode method for the quantification of the impurity was validated according to SANCO3030/99 Rev.4 and the ECHA-14-G-10-EN and the following parameters were investigated:

- Specificity and Selectivity;
- Linearity;
- Accuracy;
- Precision and Repeatability;
- LOQ

The limit content for this impurity is provided as 0.1 mg/L (100 ppb), while the value expected from the preparation of the batches was 10-50 µg/L (ppb).

For this reason, the study was set trying a range able to contain these values (40 ppb – 200 ppb).

The method for the (-)-Nicotine determination was implemented with an GC/MS - SCAN analysis in order to correctly identify the peak of analyte in the sample.

A sensitive and precise gas chromatography (GC-MS) technique was applied for determination of (-)-Nicotine. The calibration curve was found to be linear ($r = 0.9953$) over the concentration range of 0.023-0.114 µg/ml (corresponding to 0.046-0.228 µg/ml of sample). The preparations were dissolved in 2-Propanol. The samples for accuracy were prepared at 0.021, 0.053 and 0.106 µg/ml of (-)-Nicotine reference standard having known purity, representing low, middle, and high controls, respectively. Mean percentage (%) recovery \pm relative standard deviation % (RSD%) ranged from 94.78 ± 3.47 , 96.25 ± 0.03 , to 89.02 ± 0.24 . Within-day precision and instrumental repeatability were also in acceptable range: 13.8% and 1.3% respectively.

The reported method for the estimation of (-)-Nicotine proved to be specific, linear, precise, repeatable and accurate.

With this validated method were then tested 5 different production batches of the test item "VX1".

"Results" section reports the values obtained in detailed tables.

INTRODUCTION

On behalf of [REDACTED] a study aimed to validate a GC-MS method for the identification and quantification of the impurity (-)-Nicotine in the test item "VX1 validation" was performed.

Subsequently, with this validated method we tested 5 different lots of the test item production "VX1".

The study was conducted at the [REDACTED]

EXPERIMENTATION	START	END	RESEARCHER
GC/MS method for the quantification of Nicotine and five batches analysis	September 22 nd , 2016	September 24 th , 2016	[REDACTED]

BIBLIOGRAPHY

- Guidance on the Biocidal Products Regulation: Volume I: Identity/physico-chemical properties/analytical methodology - Part A: Information Requirements - Reference: ECHA-14-G-10-EN - Publ. date: November 2014.
- SANCO/3030/99 rev. 4: Technical material and preparation: Guidance for generating and reporting methods of analysis in support of pre and post registration data requirements for Annex II (part A, section 4) and Annex III (part A, section 5) of Directive 91/414.

FILING

The study program with possible amendments, raw data with possible deviations and a copy of the final report with possible revisions, will be stored in [REDACTED] archives for a period of 10 years starting from the end of the study.

At the end of the study residual sample will be kept in the fridge until December 15th, 2016, the expiry date provided by the Sponsor.

The Sponsor, drawing up of a suitable contract, may request an extension of the conservation of all or part of the documents/products for a further period or their restitution.

PROCEDURES

All procedures used during this study are recorded in the [REDACTED]

TEST ITEM IDENTIFICATION

NAME:	VX1
NATURE OF TEST SUBSTANCE:	Pesticide/Agrochemical
APPLICATION AREA:	Microbial pest control agent
STABILITY:	6 months
STORAGE:	Freezer (-18°C) and protected from light
SAMPLE DISPOSAL:	Not hazardous. The substances (plant extract) are non-toxic, non-radioactive, non-infectious and presents no risks to human or animal health or to environment.


COMPOSITION DECLARED BY THE SPONSOR :

COMPONENTS	% (w/w)
extract of tobacco plants	N.A.
dibasic sodium phosphate dodecahydrate	3
monobasic potassium phosphate	0.08
sodium sulphite	0.2
PepMV VX1	0.001-0.005

At pH 7.7 +/- 0.5

ANALYTE

NAME:	Nicotine
IUPAC NAME:	3-[(2S)-1-methylpyrrolidin-2-yl]pyridine
CAS:	54-11-5
MOLAR MASS:	162 g mol ⁻¹
FORMULA:	C ₁₀ H ₁₄ N ₂
STRUCTURE:	[Chemical structure of Nicotine: A pyridine ring attached to a pyrrolidine ring at the 3-position. The pyrrolidine ring has a methyl group (CH ₃) and a hydrogen atom (H) attached to the nitrogen atom.]



MAIN HAZARDS:

danger

RISK CODES:

H300, H310, H400, H410

ANALYZED SAMPLE FOR VALIDATION

The sample, representative of the test item, is a frozen lightly green liquid contained in a 50 ml plastic conical centrifuge tube with a screw orange cup (send 4 tubes). The liquid consists of mild isolates of pepino mosaic virus (PepMV).

Name	VX1 VALIDATION
Batch number	F
Manufacturing date	June 15 th , 2016
Expiry date	December 15 th , 2016
Receiving	EUITVI-82155
Date	Sept 02 th , 2016
#ID	ACE-2016-00123817

ANALYZED SAMPLE FOR 5 BATCHES ANALYSIS

The samples, representative of the each production batch, is a frozen lightly green liquid contained in a 15 ml plastic conical centrifuge tube with a screw blue cup (send 1 tube for each batch). The liquid consists of mild isolates of pepino mosaic virus (PepMV).

1.

Name	VX1
Batch number	A
Manufacturing date	June 15 th , 2016
Expiry date	December 15 th , 2016
Receiving	EUITVI-82155
Date	Sept 02 th , 2016
#ID	ACE-2016-00134916

2.

Name	VX1
Batch number	B
Manufacturing date	June 15 th , 2016
Expiry date	December 15 th , 2016
Receiving	EUITVI-82155
Date	Sept 02 th , 2016
#ID	ACE-2016-00134917

3.

Name	VX1
Batch number	C
Manufacturing date	June 15 th , 2016
Expiry date	December 15 th , 2016
Receiving	EUITVI-82155
Date	Sept 02 th , 2016
#ID	ACE-2016-00134918



BioPharma
Product Testing

Test Facility
Eurofins Biolab S.r.l.
GLP Provisional Certificate
Nov 20th 2015

Report No.: S-2016-03411 AM
Version: English
Page: 9 of 39
Print date: Sept 24th, 2016

4.

Name	VX1
Batch number	D
Manufacturing date	June 15 th , 2016
Expiry date	December 15 th , 2016
Receiving	EUITVI-82155
Date	Sept 02 th , 2016
#ID	ACE-2016-00134919

5.

Name	VX1
Batch number	E
Manufacturing date	June 15 th , 2016
Expiry date	December 15 th , 2016
Receiving	EUITVI-82155
Date	Sept 02 th , 2016
#ID	ACE-2016-00134920

The test item and the information concerning its was provided by the Sponsor. All data related to the test item are under the responsibility of the Sponsor and have not been verified by the Test Facility.

REFERENCE STANDARD

The (-)-Nicotine reference standard consists of a liquid contained into an amber bottle.

Name	(-)-Nicotine Pestanal (See Annex#1)
Ref. Article	[REDACTED] Supelco)
Ref No	36733
USP Batch	SZBE205XV
Assay (% w/w)	99.1
Expire date	Sept. 02 nd , 2019

Experimentation S-2016-03411 AM: GC/MS method for the quantification of Nicotine and five batches analysis

Before starting with the validation activity a set up method, in a not-GLP session, was performed in order to find a suitable method for the quantification of the impurity (-)-Nicotine in the test item VX1 validation. The optimized method was validated according to SANCO3030/99 Rev.4 guidelines and to the Guidance on the Biocidal Products Regulation ECHA-14-G-10-EN.

INFORMATION

Nicotine is a hygroscopic, oily liquid that is readily soluble in alcohol, ether or light petroleum. It is miscible with water in its base form between 60 °C and 210 °C. As a nitrogenous base, Nicotine forms salts with acids that are usually solid and water-soluble.

Nicotine is readily volatile (vapor pressure 5.5 Pa at 25°C) and dibasic ($K_{b1} = 1 \times 10^{-6}$, $K_{b2} = 1 \times 10^{-11}$).

Nicotine is optically active, having two enantiomeric forms. The naturally occurring form of Nicotine is levorotatory with a specific rotation of $[\alpha]_D = -166.4^\circ$ ((-)-Nicotine). The dextrorotatory form, (+)-Nicotine is physiologically less active than (-)-Nicotine. (-)-Nicotine is more toxic than (+)-Nicotine. The salts of (+)-Nicotine are usually dextrorotatory. The hydrochloride and sulphate salts become optically inactive if heated in a closed vessel above 180 °C.

On exposure to ultraviolet light or various oxidizing agents, Nicotine is converted to Nicotine oxide, Nicotinic acid (vitamin B3), and Methylamine.

EXPERIMENTAL PROCEDURE - VALIDATION

TEST METHOD

Gas chromatography with mass detector (GC/MS)

Parameters under investigation

Specificity>Selectivity

As part of the validation of the method, it was necessary to confirm the identity of the compound and provide that there were no relevant interferences, as required by SANCO/3030/99 rev. 4 guidelines.

For this reason, initially the solvent, the reference standard and the test sample were injected with GC/MS – SCAN Mode.

Specificity represent the capability of the method to estimate unequivocally the analyte in presence of the other components of the final product. In order to demonstrate the specificity of the method, the following solutions were separately injected into the chromatographic system:

- 1) Blank (Solvent)
- 2) (-)-Nicotine reference standard
- 3) 'VX1 validation' test sample

The specificity of the method was confirmed during the method validation.

A GC/MS-SCAN Mode confirmatory technique was used to demonstrate the method selectivity.

Linearity

Linearity refers to the ability of a detection system to produce an acceptable correlation between the instrumental response and the concentration of the analyte in the sample. The linearity of the method was assessed on the standard solutions at 40% (LOQ), 60%, 80%, 100% and 200% respectability of the limit concentration, equal to 100 µg/L. The concentrations of analyte was plotted against area. The linear regression coefficient, the slope and the intercept of the line fitting the data were calculated, along with the confidence interval at 95% for the intercept.

Accuracy

Accuracy of an analytical procedure expresses the closeness of agreement between the value which is accepted as a conventional true value and the value found with the method applied.

Two reconstituted samples (two for each level) were prepared at the concentrations corresponding to 40% (LOQ), 100% and 200% of theoretical value.

Recovery was calculated for each level and it was also determined the confidence interval of the global recovery. The accuracy was reported as mean recovery \pm relative standard deviation.

Precision

Precision of an analytical procedure refers to the closeness of agreement between mutually independent test results obtained with the same method on identical test item in the same laboratory by the same operator using the same equipment, within short intervals of time. Precision was obtained performing the assay determination of 6 samples and is expressed as RSD% of the test results.

Repeatability

Repeatability expresses the closeness of agreement (degree of scatter) between a series of measurements obtained from multiple sampling of the same homogeneous sample under prescribed conditions. Repeatability was obtained injecting a sample 6 times and is expressed as RSD% of the test results.

LOQ

The limit of quantitation (LOQ) is set at least 10 times above the blank value (expressed as Signal to Noise S/N), thus presenting a greater probability that a value at the LOQ is "real" and not just a random fluctuation of the blank reading. The LOQ is defined as the concentration at which all acceptance criteria indicated in table 1 of this study are met. The LOQ is the lowest validated level.

Acceptability criteria

The acceptability criteria for all the above-described parameters, according to SANCO/3030/99 rev. 4 and ECHA-14-G-10-EN guidelines, are summarized in the following table:

PARAMETERS	Measurement Unit	Acceptability criterion
SPECIFICITY	Chromatograms verification	No peak of blank or interferes with that of each analyte. Any interference $< 3\%$ can be neglected.
LINEARITY	R	>0.99
	Confidence limits at 95% of intercept	contain the zero
ACCURACY	% Recovery	75-125% (% impurity 'nominal' < 0.1)
	Confidence range:	$\mu = X_{average} \pm t (s/n^{1/2}) \rightarrow$ contain the 100%
PRECISION	RSD% (*)	Horwitz : $RSD_R \% = 2^{[1 - 0.5 \log(C)]}$; C is the analyte concentration as a fraction. In this case C = 0.00000010 for Nicotine (100 ppb). $RSD \% \leq Horwitz \times 0.67 = 15.2$
LOQ	% of analyte	The peak of Nicotine in the lowest solution will be a S/N ≥ 10 . This concentration falls in the linearity and accuracy range.

(*) = The acceptability of the value of RSD%, resulting from the precision, is based upon the Horwitz equation, an exponential correlation between the relative standard deviation (RSD_R) and the

concentration (C) of the analyte expressed as fraction, regardless of the analyte nature, of the matrix and of the method of measurement employed:

Horwitz equation: % RSD_R = 2^(1-0.5logC) = 22.63, considering C equal to 0.00000010

The modify Horwitz values, shown below, is used as reference in accordance with the indications of the SANCO guidelines.

$$\% RSD_r = \% RSD_R \times 0.67 = 15.2$$

Precision: RSD% ≤ RSD_r%

Analytical sequence

The analytical sequences for the quantification of the active ingredient in the precision and accuracy tests were characterized by:

- injections of two standard solutions (STD1 and STD2) containing analyte at 100% in order to evaluate the suitability of reference standard solutions to be used for the quantification;
- 6 injections of one standard solution (STD1) in order to evaluate the repeatability of the chromatographic system;
- injections of one standard solution containing analyte at 100% at the end of the analytical sequence (STD1check) in order to evaluate the suitability of complete analytical sequence;
- injections of six different sample solutions (S1, S2...S6) in order to evaluate the precision of the method.
- 6 injections of one sample solution (S1) in order to evaluate the repeatability of the method in presence of the matrix.
- injections of two reconstituted sample solution (REC SAMPLE1 and REC SAMPLE2) for each accuracy level (50%, 100% and 150%) in order to evaluate the suitability of recovery values.

Analytical acceptability criteria

The analytical acceptability criteria, according to instrumental characteristic, are summarized in the following table:


Parameters	Measurement Unit	Acceptability criteria
System suitability of standard solutions	% agreement between response factors (Fr) of standard solutions	$\% \text{ agreement STD1 vs STD2} = [\text{Ass} (F_{\text{STD1}} - F_{\text{STD2}}) / \text{average} (F_{\text{STD1}} \text{ and } F_{\text{STD2}})] \times 100 \leq 5\%$ $\% \text{ agreement STD1 vs STD1check} = [\text{Ass} (F_{\text{STD1}} - F_{\text{STD1check}}) / \text{average} (F_{\text{STD1}} \text{ and } F_{\text{STD1check}})] \times 100 \leq 5\%$
System suitability of enriched solutions	% Agreement between Recovery % values	$\% \text{ agreement REC SAMPLE 1 vs REC SAMPLE 2} = [\text{Ass} (\% \text{recovery recs1} - \% \text{recovery recs2}) / \text{average} (\% \text{recovery recs1 and } \% \text{recovery recs2})] \times 100 \leq 5\%$

The analytical sessions were considered acceptable because the criteria of the system suitability test, reported in the table above, proved to be satisfied.

Equipment

- Gas chromatograph HP 7820 A split/splitless injector, provided with MSD 5977E detector
- Capillary column Restek, Rtx®-5 Amine (30 m, 0.25 mm ID, 1.0 µm film thickness) – cat. 12353 – Serial 993086
- Standard laboratory equipment

The GC instrument used is qualified every year according internal procedure.



Reagent

- 2-propanol - IPA (Sigma-Aldrich, code 190764, lot n. BCBR0008V)

All standard and reagents are of high purity analytical grade. The validity of each product was checked before starting the analyses.

Working solutions preparation

Specificity

Blank

2-Propanol.

(-)-Nicotine reference standard - mother solution (conc. ~ 2 mg/ml)

About 100 mg of (-)-Nicotine reference standard were quantitatively weighed into a 50 ml volumetric flask and brought to volume with 2-Propanol (Solution A).

(-)-Nicotine reference standard - intermediate solution (conc. ~ 0.01 mg/ml)

100 µl of (-)-Nicotine reference standard mother solution were quantitatively transferred into a 20 ml volumetric flask and brought to volume with 2-Propanol (Solution B).

(-)-Nicotine reference standard - LOQ solution (conc. ~ 0.02 µg/ml)

40 µl of (-)-Nicotine reference standard (solution B) were quantitatively transferred into a 20 ml volumetric flask and brought to volume with 2-Propanol.

Test sample

About 1.0 ml of test sample was quantitatively weighed into a 2.0 ml volumetric flask and brought to volume with 2-Propanol.

Linearity

(-)-Nicotine reference standard - mother solution (conc. ~ 2 mg/ml)

About 100 mg of (-)-Nicotine reference standard were quantitatively weighed into a 50 ml volumetric flask and brought to volume with 2-Propanol (Solution A).

(-)-Nicotine reference standard - intermediate solution (conc. ~ 0.01 mg/ml)

100 µl of (-)-Nicotine reference standard mother solution were quantitatively transferred into a 20 ml volumetric flask and brought to volume with 2-Propanol (Solution B).

Calibration:

Level 30% - 0.03 µg/ml of sample	30 µl of Nicotine mother solution B	to 20 ml with 2-Propanol
Level 60% - 0.06 µg/ml of sample	60 µl of Nicotine mother solution B	to 20 ml with 2-Propanol
Level 80% - 0.08 µg/ml of sample	80 µl of Nicotine mother solution B	to 20 ml with 2-Propanol
Level 100% - 0.1 µg/ml of sample	100 µl of Nicotine mother solution B	to 20 ml with 2-Propanol
Level 200% - 0.2 µg/ml of sample	200 µl of Nicotine mother solution B	to 20 ml with 2-Propanol

Precision/Repeatability

Sample preparation (six preparations)

About 1.0 ml of test sample was quantitatively weighed into a 2.0 ml volumetric flask and brought to volume with 2-Propanol.

The sample must be thawed slowly in refrigerator (5 ± 3 °C) and prepared just before the injections.

Precision was obtained performing the assay determination of these 6 samples and is expressed as RSD% of the test results.

Repeatability was obtained injecting the first sample 6 times and is expressed as RSD% of the test results.

Accuracy

(-)-Nicotine reference standard - mother solution (conc. ~ 2 mg/ml)

About 100 mg of (-)-Nicotine reference standard were quantitatively weighed into a 50 ml volumetric flask and brought to volume with 2-Propanol (Solution A).

(-)-Nicotine reference standard - intermediate solution (conc. ~ 0.01 mg/ml)

100 µl of (-)-Nicotine reference standard mother solution were quantitatively transferred into a 20 ml volumetric flask and brought to volume with 2-Propanol (Solution B).

Calibration:

STD - Level 1 (40%)	40 µl of Nicotine mother solution B	to 20 ml with 2-Propanol
STD - Level 2 (60%)	60 µl of Nicotine mother solution B	to 20 ml with 2-Propanol
STD - Level 3 (80%)	80 µl of Nicotine mother solution B	to 20 ml with 2-Propanol
STD - Level 4 (100%)	100 µl of Nicotine mother solution B	to 20 ml with 2-Propanol
STD - Level 5 (200%)	200 µl of Nicotine mother solution B	to 20 ml with 2-Propanol

Enriched sample 40% (two preparations)

1.0 ml of test sample was quantitatively weighed into a 2.0 ml volumetric flask. 1.0 ml of Reference standard - Level 3 was quantitatively transferred and the solution was brought to volume with 2-Propanol.

Enriched sample 100% (two preparations)

1.0 ml of test sample was quantitatively weighed into a 2.0 ml volumetric flask. 1.0 ml of Reference standard - Level 5 was quantitatively transferred and the solution was brought to volume with 2-Propanol.

Enriched sample 200% (two preparations)

1.0 ml of test sample was quantitatively weighed into a 2.0 ml volumetric flask. 20 µl of Reference standard - Solution B were quantitatively transferred and the solution was brought to volume with 2-Propanol.

LOQ

The LOQ is the analyte concentration with S/N ratio of at least 10 and was verified on the working standard solutions at 0.02 µg/ml.

To conduct excellently this delicate parameter, in addition to verify that this level was in the linearity range and met the accuracy in presence of matrix, we were carried out three other preparations starting from different mothers to measure the reproducibility of its response.

(-)-Nicotine reference standard - mother solution (conc. ~ 2 mg/ml) – Three preparations

About 100 mg of (-)-Nicotine reference standard were quantitatively weighed into a 50 ml volumetric flask and brought to volume with 2-Propanol (Solution A).

(-)-Nicotine reference standard - intermediate solution (conc. ~ 0.01 mg/ml) – Three preparations

100 µl of (-)-Nicotine reference standard mother solution were quantitatively transferred into a 20 ml volumetric flask and brought to volume with 2-Propanol (Solution B).

(-)Nicotine reference standard - LOQ solution (conc. ~ 0.02 µg/ml) – Three preparations
40 µl of (-)-Nicotine reference standard - solution B were quantitatively transferred into a 20 ml volumetric flask and brought to volume with 2-Propanol.

Five batches analysis

As required by Guidance on Regulation (EU) No 528/2012 guidelines, it was necessary to determine the analytical profile of five representative production batches of test item.
For this reason, five production batches were analysed, confirming the repeatability of the production process.

Sample preparation

About 1.0 ml of test sample was quantitatively weighed into a 2.0 ml volumetric flask and brought to volume with 2-Propanol.
The sample must be thawed slowly in refrigerator (5 ± 3 °C) and prepared just before the injections.

GC-MS METHOD – SIM MODE

According to good mass spectrometric analysis, a minimum of 3 ions (ideally with an m/z ratio of >100) must be used for identification/quantification.

Instrumentation	GC-MS
Column	Rtx®-5 Amine, 30m x 0.25mm x 1.0 µm
Detector (Aux 2)	MSD, 280°C
Source	230°C
Quadrupole	150°C
MS mode	SIM (84 m/z, 161 m/z, 162 m/z)
Threshold	100
Flow	Helium, 1.1 ml/min
GC oven program	100°C for 0 min, rate 20°C/min to 200°C for 0 min, rate 35°C/min to 300°C for 5 min.
Run Time	12.86 min
Injector temperature	290°C
Injection volume	1 µl - Split 1:2
Solvent Delay	4.00 min
RETENTION TIME	(-)Nicotine ~ 6.2 min

CALCULATIONS

The content of analyte was calculated as follows:

$$C_{ST} = \frac{W_{ST} \times T_{ST}}{D_{ST} \times 100} \times 1000 ;$$

Where:

C_{ST} = standard concentration (mg/ml)

W_{ST} = standard weight (g)

D_{ST} = dilution of analyte in working standard solution (ml)

T_{ST} = standard assay (% w/w)

A calibration line was obtained after injections of diluted reference standard solutions at the concentrations of 0.02, 0.03, 0.04, 0.05 and 0.1 µg/ml (respectively 0.04, 0.06, 0.08, 0.1 and 0.2 µg/ml of sample) :

$$Y = aX + b$$

$$(-)-NICOTINE (\mu\text{g/ml}) = [(Y - b)/a] \times D_c$$

Where:

$Y = A_{ST}$ = Area of the analyte in working standard solution

$X = C_{ST}$ = Conc. (µg/ml) of the analyte in working standard solution

A_{ST} = analyte area in the working standard solution (pA*s)

a = slope

b = intercept

$$(-)-NICOTINE (\% \mu\text{g}/\mu\text{g}) = \frac{[(Y - b)/a] \times D_c}{W_c} \times 100$$

Where:

D_c = sample dilution (ml)

W_c = sample weight (µg)

The quantifications are made on the 84 m/z fragment, being the most abundant. The other two fragmentation ions are measured as verification of the abundances.

EXPERIMENTAL PROCEDURE - GC/MS CONFIRMATORY TECHNIQUE

TEST METHOD – GC/MS

Gas chromatography with mass detector (GC/MS)

Equipment

- Gas chromatograph HP 7820 A split/splitless injector, provided with MSD 5977E detector
- Capillary column Restek, Rtx®-5 Amine (30 m, 0.25 mm ID, 1.0 µm film thickness) – cat. 12353 – Serial 993086
- Standard laboratory equipment

The GC instrument used is qualified every year according internal procedure.

Analysis

The method was implemented with an analysis GC/MS in order to correctly identify the peaks of the analyte in the sample and in the raw material.

For this reason the standard solution (Solution B) and the test sample preparations were injected with GC/MS technique - SCAN mode to confirm the analyte identity, while the analysis was carried out and validate in SIM mode (after selection of the characteristic fragments).

Working solutions preparation

Blank

2-Propanol.

(-)-Nicotine reference standard - mother solution (conc. ~ 2 mg/ml)

About 50 mg of (-)-Nicotine reference standard were quantitatively weighed into a 50 ml volumetric flask and brought to volume with 2-Propanol (Solution A).

(-)-Nicotine reference standard - intermediate solution (conc. ~ 0.01 mg/ml)

100 µl of (-)-Nicotine reference standard mother solution were quantitatively transferred into a 20 ml volumetric flask and brought to volume with 2-Propanol (Solution B).

Test sample

About 1.0 ml of test sample was quantitatively weighed into a 2.0 ml volumetric flask and brought to volume with 2-Propanol.

These solutions were injected using the following GC/MS method.

Instrumentation	GC-MS
Column	Rtx®-5 Amine, 30m x 0.25mm x 1.0 µm
Detector (Aux 2)	MSD, 280°C
Source	230°C
Quadrupole	150°C
MS mode	SCAN (30-200 amu)
Threshold	100
Flow	Helium, 1.1 ml/min
GC oven program	100°C for 0 min, rate 20°C/min to 200°C for 0 min, rate 35°C/min to 300°C for 5 min.
Injector temperature	290°C

Injection volume	1 μ l - Split 1:2
Solvent delay	4.00 min


The following parameters were evaluated:

- 1) characteristic MS ions: characteristic MS ions were extrapolated by mass spectra.
- 2) identification: identity of each peak were identified by mass spectra analyses and interpretation of the MS spectra performed by using the NIST/EPA/NIH (version-11.0).
- 3) Probability %: the correspondence between the mass spectra obtained by analysis and those contained in the NIST/EPA/NIH Library version-11.0 is expressed through the parameter "probability %." This value indicates how the unknown substance is correctly identified from the reference library. Values greater than 90% indicate a good correlation, while values below 50% indicate that there is a substantial difference between the compound analysed and the reference library. Differences of \pm 5% in the values of probability are not considered significant.

According to the conditions previously described, the peaks belonging to Nicotine was clearly visible. The recognition by the GC/MS library has allowed to characterize this principal peak present in the solutions. This peak were under mentioned. In the range time (4.0 min) the solvent (2-propanol) fell and to avoid mass source damage the solvent delay was set.

For this reason, in this zone, other possible compounds were not recognized.

In the following table, identifications of the investigated compound is summarized.



Group	Retention Time (min)	Selected Ions ^a			
		1	2	3	4
Nicotine	1.7	84 (100)	133 (22) ^b	161 (16)	162 (M) ^c (18)

^a Dwell times should be adjusted to produce a cycle time of about 4 scans/sec

^b Percent relative abundance with respect to the ion of highest abundance.

^c M is the Molecular Ion.

RESULTS

VALIDATION METHOD PARAMETERS

Selectivity

The method was implemented with an analysis GC/MS in order to correctly identify the peak of the analyte Nicotine in the reference standard and in the sample. The GC/MS-SCAN Mode has proved itself as confirmatory technique.

Figure 1: Blank mass chromatogram




Figure 2: Test sample mass chromatogram







Figure 3: (-)-Nicotine mass chromatogram in the reference standard (0.01 mg/ml)




Library Searched : D:\MassHunter\Library\NIST11.L
Quality : 96
ID : Pyridine, 3-(1-methyl-2-pyrrolidinyl)-, (S)-



Specificity

The method proved to be specific; in fact it has been verified that the test solution do not interfere with the peak of Nicotine.

Figure 1: Blank chromatogram (Solvent)



Since mass chromatogram was normalized on the major peak, a narrower time range (4.0 – 9.0 min) was selected to better highlight the peak of interest ($rt \approx 6.2$ min).

Figure 2: Sample chromatogram




Figure 3: Nicotine chromatogram (Reference standard 0.02 µg/ml – LOQ level)




Figure 4: Nicotine chromatogram (Reference standard 0.05 µg/ml – 100% level)






Figure 5: Enriched sample chromatogram (LOQ addition)



Linearity

Method Linearity was tested on the Nicotine reference standard on 5 different concentration levels from 30% to 200% of the theoretical amount of analyte in the sample.
 All acceptance criteria ($R > 0.99$ and/or the confidence interval at 95% for the intercept contains zero) were satisfied (see Annex#2: excel sheet).



(-)-Nicotine Linearity	
R	0.9931
Slope	1168127.2195
Intercept	8947.6376
Confidence interval at 95% for the intercept	[-0.0088;0.0088]

Accuracy

Method accuracy was tested on spiked samples prepared at three concentration levels for analyte (40% (LOQ), 100% and 200% of the theoretical amount of analyte in the sample). Two fortified samples (two for each level) were prepared (see Annex#4: excel sheet).

(-)-Nicotine				
	Recovery %	Average Recovery %	95% confidence interval lower limit	95% confidence interval upper limit
Reconstituted sample at 40%	99.02			
Reconstituted sample at 100%	102.49	102.82	98.76	106.89
Reconstituted sample at 200%	106.97			

The measured values were compared with the 'expected' value of 100% using the Student's t-test and the choice of null hypothesis was appropriate to the data set.

For the t-tests the following Equation A was used:

Equation A:

$$t_{\text{cal}} = \frac{\bar{x} - \mu}{s / \sqrt{n}}$$

where

\bar{x} = mean of test results of a sample

μ = "true" or reference value

s = standard deviation of test results

n = number of test results of the sample.

To compare the mean of a data set with a reference value normally the "two-sided t-table of critical values" is used ($t_{\text{cal}} \leq t_{\text{tab}}$). The applicable number of degrees of freedom here is:

$$df = n - 1$$

The value for t calculated with Equation A did not exceed the critical value in the table, therefore the data were taken to belong to the same population: there is no difference and the "null hypothesis" is accepted (with the applicable probability of 95%).

The acceptance criterion (Recovery active ingredient = 95%-105%, % and/or the confidence interval at 95% for the recovery contains 100%) proved to be satisfied.

Precision

Method precision was proved preparing 6 different samples in the same analytic session. The RSD% of the percentage assay (% w/w) was calculated for the analyte on the test sample. Results are reported below (see Annex#3: excel sheet). The experimental RSD% respects the acceptance criteria.

Analyte	Assay (% w/w)	Acceptance criteria	RSD%
Nicotine content	0.0000072	RSD% NICOTINE \leq 15.2% ^(*)	1.3

(*) = The Horwitz equation is an exponential correlation between the relative standard deviation (RSD_R) and the concentration (C) of the analyte expressed as fraction, regardless of the analyte nature, of the matrix and of the method of measurement employed:

$$\% \text{ RSD}_R = 2^{(1-0.5\log C)}$$

Repeatability

Repeatability was obtained injecting a sample 6 times and is expressed as RSD% of the test results.

Results are reported below.

The method proved to be repeatable (see Annex#3: excel sheet).

Analyte in the preparation S1	Assay (% w/w)	RSD%
Nicotine content	0.0000060	1.7

As you can see the sample undergoes degradation over time. Repeatability is fact was proven by injecting the first vial of precision, put on the GC sampler.

LOQ

The LOQ is defined as the concentration at which all acceptance criteria indicated in table "Acceptability criteria" (page 11 of this report) of this study are met. The LOQ is the analyte concentration at which the S/N ratio is at least 10 and corresponds to the lowest validated level

The method proved that the concentration corresponding to 0.02 µg/ml (or 0.04 µg/ml on the sample) had a signal to noise ≈ 30 and fall in the linearity and accuracy range (see Annex#6: excel sheets).

ANOMALY TEST (DIXON, GRUBBS):

The precision data set were subjected to statistical analysis to verify the presence of values statistically abnormal.

On results, originated from the same sample or from similar samples, it is necessary to verify the adequacy of the data to be statistically processed by appropriate tests that allow us to discard outliers.

In this regard we have carried out the Dixon and Grubbs tests, formulating the null hypothesis H0 as: the data is not to be excluded from the calculations.

For none of the tests the null hypothesis was rejected (see Annex#7: excel sheets).

Statistical consideration

The measured precision is within the recommended values, given by Horwitz modification values.

It is then possible to derive the reproducibility standard deviation, σ_R , from the approximate form of the Horwitz equation.

$\sigma_R = 0.02 * C^{0.8495}$ which as you can see directly puts in relation σ_R with the analyte concentration.

Next, multiplying σ_R for the coverage factor, you get the expanded uncertainty.

Before using the reproducibility standard deviation for the calculation of the expanded uncertainty, it is necessary to verify that its close repeatability standard deviation (S_r) is compatible with σ_R obtained from the Horwitz equation.

It have to check the condition $1/2 \sigma_R \leq S_r \leq 2/3 \sigma_R$.

You can have a better repeatability standard deviation, occurring $1/2 \sigma_R > S_r$.

In our method there is the condition in which our values are below the lower limit ($\sigma_R \times 0.5$).

% w/w	Nicotine
Xmedium	0.0000072
$U_e = K * \sigma_R$	Non applicable

However, since our analyte content is affected by errors, in excess or defect, we prefer to give an estimate interval that expresses this error in the equation below (obtained from 6 preparations - see Annex#5: excel sheet):

$$X = X_{\text{medium}} \pm t^* S_r / RADQ(n)$$

% w/w	Nicotine
Xmedium	0.0000072
± t*Sr/RADQ(n)	± 0.0000001

DETERMINATION OF (-)-NICOTINE CONTENT IN FIVE BATCHES OF THE TEST ITEM "VX1"

Sample ID	Batch	(-)-Nicotine (*)
ACE-2016-00134916	A	0,069 µg/ml
ACE-2016-00134917	B	0,073 µg/ml
ACE-2016-00134918	C	0,072 µg/ml
ACE-2016-00134919	D	0,080 µg/ml
ACE-2016-00134920	E	0,082 µg/ml
	Average (µg/ml)	0,075
	SD (µg/ml)	0,005
	RSD%	7.2

(*) Nicotine LOQ = 0.04 µg/ml = 0,000004 % w/w, according to Validated method S-2016-03411AM.
 See Annex#4 for individual data and calculations.

DEVIATION

No deviation has been recorded from study program.

CONCLUSIONS

The method described in this study proved to be specific, linear, precise and repeatable and was successfully validated.

The (-)-Nicotine average content in five production batches of test item was 0.075 µg/ml with a deviation standard equal to 0.005 µg/ml.

ANNEXES

ANNEX	TITLE
N.1	NICOTINE - REFERENCE STANDARD CoA
N.2	LINEARITY – EXCEL SHEET
N.3	PRECISION-REPEATABILITY – EXCEL SHEET
N.4	ACCURACY AND 5 BATCHES ANALYSIS – EXCEL SHEET
N.5	STATISTIC – EXCEL SHEET
N.6	LOQ – EXCEL SHEET
N.7	ANOMALY TEST – EXCEL SHEET

ANNEX#1: NICOTINE - REFERENCE STANDARD CoA

CERTIFICATE OF ANALYSIS

Sigma-Aldrich Laborchemikalien GmbH D-30918 Seelze
Telefon: +49 5137 8238-150

Seelze, 17.09.2014/541285/14/12728

Order-No.:
Customer-No.:

Order-Code:

Quantity:

Production Date: 02.Sep.2014
Expiry Date: 02.Sep.2019

Article/Product: 36733

Batch : SZBE205XV

(-) -Nicotine PESTANAL®

Reference Material (RM)

1. General Information

Formula: C₁₀H₁₄N₂
CAS-No.: [54-11-5]
Usage : Insecticide

Molar mass: 162.23 g/Mole
Recomm. storage temp.: roomtemp.

The estimated uncertainty of a single measurement of the assay can be expected to be 0.5 % relative (confidence level = 95%, n= 6) whereby the assay measurements are calculated by 100% minus found impurities.

2. Batch Analysis

Identity (NMR)

complying

Assay (GC)

99.1 area %

Refractive index (n 20/D)

1.5278

Date of Analysis

17.Sep.2014

3. Advice and Remarks


- The expiry date is based on the current knowledge and holds only for proper storage conditions in the originally closed flasks/ packages.
- Whenever the container is opened for removal of aliquot portions of the substance, the person handling the substance must assure, that the integrity of the substance is maintained and proper records of all its handlings are kept. Special care has to be taken to avoid any contamination or adulteration of the substance.
- We herewith confirm that the delivery is effected according to the technical delivery conditions agreed.
- Particular properties of the products or the suitability for a particular area of application are not assured.
- We guarantee a proper quality within our General Conditions of Sales.

GLC-Method

Analytical Department

Article : (-)-Nicotine
Article-No : 36733
Batch : SZBE205XV
Column : SP-1701, 30m, 0,32mm i.D., 1.0µm Film
Inj.-Temp. : 280°C
Det.-Temp. : 280°C - FID
Oven-Temp. : 150°C to 250°C (10°C/min) hold 20min
Split : 1:100
Flow : 1ml He/min
Inj.v. : 0,2µl
Evaluation : uncorrected
Operator : Schulz

C H R O M A T O G R A M M



Area Percent Report

#	Meas.	Re	Height	Area	Area %
1		8.60	8886.7	30157.3	99.16
2		8.96	1.0	5.2	0.02
3		9.19	7.1	27.1	0.09
4		9.92	3.1	8.2	0.03
5		10.00	35.6	99.0	0.33
6		10.25	1.1	5.0	0.02
7		10.61	2.5	13.7	0.04
8		13.31	5.0	19.2	0.06
9		16.56	11.1	61.8	0.20
10		19.49	2.3	16.1	0.05


2/3
23/09/16 UF

Report No.: S-2016-03411 AM
Version: English
Page: 29 of 39
Print date: Sept 24th, 2016

(-) -Nicotin
PSTANAL
#36733 Ch. :SZ3E2Q5XY

Sample Name:
14_18728
Data Collected on:
Sesize-NMR vnmr346Q
Archive directory:
/home/vnmr1/vnmrsys/data
Sample directory:
14_18728_20140915
FidFile: 14_18728_PROTON_01
Pulse Sequence: PROTON (s2pnl)
Solvent: cdd13
Data collected on: Sep 15 2014


Temp. 25.0 C / 299.1 K
Sample #2, Operator: vnmr1
Relax. delay 5.000 sec
Pulse 45.0 degrees
Acq. time 2.261 sec
Width 7183.9 Hz
16 repetitions
OBSERVE RL, 399.8818069 MHz
DATA PROCESSING
FT size 32768
Total time 1 min 56 sec



3/3
23/09/16 VP

Report No.: S-2016-03411 AM
 Version: English
 Page: 30 of 39
 Print date: September 24th, 2016

ANNEX#2 : LINEARITY - EXCEL SHEET



Report No.: S-2016-03411 AM
 Version: English
 Page: 31 of 39
 Print date: Sept 24th, 2016

ANNEX#3 : PRECISION-REPEATABILITY - EXCEL SHEET

STUDY NUMBER		S-2016-03411 AM			
DATE	24-Sep-16	PARAMETER TO VALIDATE	PRECISION AND REPEATABILITY		
PRODUCT	Vx1				
STANDARD MOTHER SOLUTIONS:					
Reference	Weight (g)	Dilution (ml)	Cone. 1 (mg/ml)		
Nicotine	0.1051	99.1	2.083		
Nicotine	0.1160	50	2.299		
Reference	Sampling (ml)	Dilution (ml)	Cone. 2 (mg/ml)		
STD 1	0.1	20.0	0.010		
STD 2	0.1	20.0	0.011		
STANDARD JOB SOLUTIONS:					
STD Level	Sampling (l)	STD Concentration (ug/ml)	STD Area (STD/Area) (average)		
Level 40%	30.0	0.016	23411.5		
Level 60%	60.0	0.031	43678		
Level 80%	80.0	0.042	62267		
Level 100%	100.0	0.052	59538		
Level 200%	200.0	0.104	72889		
Slope	1188127.2195		73232		
Intercept (R)	8947.6376		73060.5		
Correlation (R ²)	0.9966		7.14E-07		
R ²	0.9931		7.11E-07		
PRECISION					
Sample preparations	Weight (g)	Sample Dilution Volume (ml)	Nicotine conc. (ug/ml)		
1	1.0081	(R= 6.20 min, m/z 84) 50938	0.072		
2	1.0110	51578	0.073		
3	1.0100	52157	0.074		
4	1.0076	51063	0.072		
5	1.0088	51194	0.072		
6	1.0087	52314	0.074		
		51541	0.073		
		1.1	0.0000072		
Standard Solution	Sampling (l)	Volumetric flask (ml)	Conc. (mg/ml)		
STD 2	100.0	20.0	0.057		
STANDARD SOLUTION 2		Area (Conc/Area) (Average)	Fr (Averag)		
STD 2		81452	7.06E-07		
		78663	7.32E-07		
		Assessment -in (%)	Assessment -in (%)		
STD1 - STD2		3.61	0.82		

	Report No.: S-2016-03411 AM
Version:	English
Page:	32 of 39
Print date:	Sept 24 th , 2016

REPEATABILITY		Sample preparations		Sample Dilution Volume (ml)		Nicotine conc. (µg/ml)		Nicotine %	
1	Weight (g)	Area (R= 6.20 min, m/z 84)							
2		44273						0.00000595	
3	1.0163	43519						0.00000582	
4		44944						0.00000606	
5		44061		2.0				0.00000582	
6		45188						0.00000611	
		44697						0.00000602	
		44447						0.0000060	Average
		1.4						0.0000060	RSD%
								1.7	

STANDARD SOLUTION 1 check		Sampling (µ)		Volumetric flask (ml)		Conc. (ng/ml)		Fr (Average) (%)		Assessment-inj (%)		Assessment	
Standard Solution												STD1 - STD1check	
STD 1 check		100.0		20.0		0.052		7.05E-07	7.18E-07	7.30E-07	3.53	0.69	

Report No.: S-2016-03411 AM
 Version: English
 Page: 33 of 39
 Print date: Sept 24th, 2016

ANNEX#4 : ACCURACY AND 5 BATCHES ANALYSIS - EXCEL SHEET

STUDY NUMBER		S-2016-03411 AM	
DATE		24-Sep-16	
PARAMETER TO VALIDATE		ACURACY	
PRODUCT		Vx-1	
STANDARD MOTHER SOLUTIONS:			
Reference	Weigh (g)	Dilution (ml)	Conc. 1 (mg/ml)
Nicotine	0.1040	50	2.081
Nicotine	0.1160	50	2.299
STANDARD JOB SOLUTIONS:			
Reference	Sampling (ml)	Dilution (ml)	Conc. 2 (mg/ml)
STD 1	0.1	20.0	0.010
STD 2	0.1	20.0	0.011
STD Level	Sampling (ml)	Dilution (ml)	STD Concentration (μg/ml)
Level 40%	40.0	20.0	0.021
Level 60%	60.0	20.0	0.031
Level 80%	80.0	20.0	0.041
Level 100%	100.0	20.0	0.052
Level 200%	200.0	20.0	0.103
Slope	968562.6143		
Intercept	5941.0515		
Correlation (R)	0.9926		
R ²	0.9856		
Value to subtract (Area)	44447		
ENRICHED SOLUTION			
Sample preparations	Weight (g)	Theoretical conc. (μg/ml)	Sample Dilution Volume (ml)
LOQ-1	1.0082	0.021	(Re= 6.20 min, m/z 84)
LOQ-2	1.0083	0.021	64398 19851 0.021
100%-1	1.0081	0.052	65313 20866 0.021
100%-2	1.0067	0.052	105199 60752 0.052
200%-1	1.0083	0.103	105034 61487 0.053
200%-2	1.0082	0.103	155629 111482 0.111
			154550 110103 0.110
Confidence interval:	t _{0.95} (n=11; 95%)	RADQN	Recovery (%)
t = Xaverrate ± t(s/n) ^a	2.571	2.45	98.76
t ^b Tests for bias:	X average - H	s ^c t calc	Assessment % (between preparations)
	2.82	9.48	4.48
GLOBAL RECOVERY			
% w/w average	102.62		SYSTEM SUITABILITY (%)
Standard deviation	3.87		Agreement % of STD1-STD2
Global RSD%	3.77		Agreement % of STD1-Dcheck
			Agreement % between recovery level
			VALID
			VALID
			VALID



Report No.: S-2016-03411 AM
 Version: English
 Page: 34 of 39
 Print date: Sept 24th, 2016

FIVE BATCHES ANALYSIS

Sample Preparations	Weight (g)	Area ($R=6.20 \text{ min, m/z 84}$)	Sample Dilution Volume (ml)	Analyte Assay conc. (ug/ml)	Analyte Assay %
VC-1 lot A	1.0084	364693		0.059	0.000069
VC-1 lot B	1.0087	41167		0.073	0.000072
VC-1 lot C	1.0074	40759	2.0	0.072	0.000072
VC-1 lot D	1.0085	44395		0.060	0.000079
VC-1 lot E	1.0084	456658		0.062	0.000081
				Average	0.000076
				SD	0.000005
				RSD%	7.1

STANDARD SOLUTION 2

Standard Solution	Sampling (μl)	Volumetric flask (ml)	Cone.	Area ($R=6.20 \text{ min, m/z 84}$)	Fr (Conc/Area)	Assessment ini (%)	Assessment std (%)
STD 2	100.0	20.0	0.057	69494	8.21E-07	6.98	0.83
				64808	8.81E-07		

	Report No.: S-2016-03411 AM
Version:	English
Page:	35 of 39
Print date:	Sept 24 th , 2016

ANNEX#5 : STATISTIC - EXCEL SHEET

STUDY NUMBER		S-2016-03411 AM	
DATE		23 September 2016	
PARAMETER TO VALIDATE		HORWITZ	
PRODUCT		VX-1	
Sample ID		Unit	Set data
ACE-2016-00123817		% w/w	0,0000071 0,0000072 0,0000073 0,0000072 0,0000072 0,0000074
t student (two sided):		n-1	95% 99%
		2	4,303 9,925
		3	3,182 5,841
		4	2,776 4,604
		5	2,571 4,032
		6	2,447 3,707
		7	2,365 3,499
		8	2,306 3,356
		9	2,262 3,250
		10	2,228 3,169
		11	2,201 3,106
number of tests		6	
t value (95%, n-1)		2,571	% w/w
Media		0,0000072	% w/w
Standard deviation (Sr)		9,47E-08	% w/w
Variance (S^2)		8,98E-15	% w/w
Limit of repeatability r		0,0000003	% w/w
RSD%			
theoretical value			
Definition of C (mass/mass)			
Horwitz equation (% RSD _R):			
$\sigma_R (\%)$:			
Horwitz corrected (% RSDr) :			
RESULTS:		X = Xmedium ± U _E X = Xmedium ± t*Sr/RADQ(n)	0,0000072 ± Not applicable 0,0000072 ± 0,0000001
U _E = K * σR =			
Acceptedability of the value of RSDr%:		RSD% ≤ %RSDr	VALID
Significance repeatability 1/2 σR ≤ Sr ≤ 2/3 σR:		1/2 σR ≤ Sr	≤ 2/3 σR
Expanded uncertainty:		0,00000	0,0000

Report No:	S-2016-03411 AM
Version:	English
Page:	36 of 39
Print date:	Sept 24 th , 2016

ANNEX#6 : LOQ - EXCEL SHEET

STUDY NUMBER		S-2016-03411 AM	
DATE		24-Sep-16	
PARAMETER TO VALIDATE		LOQ	
PRODUCT		VX-1	
STANDARD MOTHER SOLUTIONS:			
Reference	Weigh (g)	Assay (%)	Dilution (ml)
Nicotine	0,1160	99,1	50
Nicotine	0,1080	99,1	50
Nicotine	0,1065	99,1	50
Sampling (ml)	Dilution (ml)	Conc. 1 (mg/ml)	
STD 1	0,1	20,0	0,011
STD 2	0,1	20,0	0,011
STD 3	0,1	20,0	0,011
STANDARD JOB SOLUTIONS:			
STD Level	Sampling (μ l)	STD Concentration (μ g/ml)	STD Area (R= 6,20 min, m/z 84) (average)
40% (LOQ) n.1	40,0	0,023	22528
40% (LOQ) n.2	40,0	0,021	21764
40% (LOQ) n.3	40,0	0,021	20672
		Average RSD%	
		2,1	

	Report No.:	S-2016-03411 AM
	Version:	English
	Page:	38 of 39
	Print date:	Sept 24 th , 2016

2) Grubbs Test

STUDY NUMBER	DATE	ANOMALY TEST						
		PARAMETER TO VALIDATE	PRODUCT	ANOMALY TEST - GRUBBS (UNICHIM n.179/2)				
Id. number	set	copy	Significance levels:					
			L	S	1%	5%		
ACE-2016-00123817			$G_{1,n}$ (Largest)	1,411	Correct value	4	1,496	1,481
			$G_{1,1}$ (Smallest)	1,009	Correct value	5	1,764	1,715
						6	1,973	1,887
						7	2,139	2,020
						8	2,274	2,126
			Critical values G_2	5%	1%	9	2,387	2,215
			Sk,Lk	1,887	1,973	10	2,482	2,290
						11	2,564	2,355
						12	2,636	2,412
			number of suspicious values at each end					
			n	6	1			
			s	0,000				
			Media	0,000				

NOTE: The extreme value (min/max) can be considered:

- "Correct", and then accepted as consistent with the other data, if the criterion is lower than the significance level of 5%.
- "Dispersed", if the criterion met the significance level of 1%.
- "Anomalous", if the criterion does not meet any significance level.

